



25470



25468



24986



25469



25471

Injector & Autosampler Supplies for Shimadzu HPLC Systems

| Description | Model # | Similar to Shimadzu part # | qty. | cat.# |
|---------------------|----------------------------|----------------------------|------|-------|
| Needle Seal, Vespel | LC-2010, SIL-10ADvp, 10AXL | 228-33355-04 | ea. | 25468 |
| Rotor Seal | SIL-10ADvp | 228-21217-97 | ea. | 24986 |
| Rotor Seal Assembly | SIL-10A, 10AXL, 10Ai | 228-21217-91 | ea. | 25469 |
| Stator Assembly | SIL-10A, 10AXL, 10Ai | 228-21220-91 | ea. | 25470 |
| Syringe, 500 µL | SIL-10A, 10AXL, 10Ai | 228-25237-04 | ea. | 25471 |

For even more options, visit www.restek.com/lcacc




25284

Deuterium Lamps for Shimadzu HPLC Systems

| Description | Model # | Similar to Shimadzu part # | qty. | cat.# |
|----------------|--|----------------------------|------|-------|
| Deuterium Lamp | SPD-10, 10A, 10AVvp, SPD-20A, SPD-20AV | 228-34016-02 | ea. | 25284 |

Protect your column and your UHPLC performance with UltraShield and UltraLine UHPLC Filters

A cost-effective way to extend the lifetime of any UHPLC column without sacrificing your UHPLC performance on any LC system.



www.restek.com/LCguard

Instrument Replacement Parts for Shimadzu HPLC Systems




Pump Supplies for Shimadzu HPLC Systems

| Description | Model # | Similar to Shimadzu part # | qty. | cat.# |
|---|--|----------------------------|------|-------|
| Plunger Seal, Polyethylene, Gold** | LC-10Ai, 10AS, 10AT, 10ATvp, LC-7A | 228-21975-00 | ea. | 25290 |
| Plunger Seal, Black* | LC-10AD, 10ADvp, LC-20AD/AB, LC-600, LC-2010 A/C HT, LC-9A | 228-35146-00 | ea. | 24980 |
| Plunger Seal, Gold** | LC-10ADvp, LC-20AD/AB | 228-32628-00 | ea. | 24981 |
| Plunger Seal, Black* | LC-10ATvp, LC-20AT, LC-HT, SIL-10ADvp, SIL-2010 A/C HT, SIL-20A/AC, SIL-HT | 228-35145-00 | ea. | 24985 |
| Plunger Rinse Seal | LC-10Ai, LC-10AS, LC-10AT, LC-10ATvp | 228-28499-00 | ea. | 25292 |
| Inlet Check Valve | LC-6A, LC-10AS | 228-12353-91 | ea. | 25287 |
| Inlet Check Valve | LC-600, LC-9A, LC-10AD | 228-18522-91, 228-33492-91 | ea. | 25295 |
| Inlet Check Valve | LC-10ADvp | 228-39093-92 | ea. | 24984 |
| Inlet Check Valve | LC-10AT, LC-10ATvp | 228-32166-91 | ea. | 26521 |
| Inlet Check Valve | LC-20AD/AB XR | 228-48249-91, 228-45704-91 | ea. | 26426 |
| Outlet Check Valve | LC-6, LC-10AS, LC-8A | 228-09054-93 | ea. | 25288 |
| Outlet Check Valve Rebuild Kit Includes: ruby ball (2), seat (2), thin washer (2), thick washer (2), retainer washer (2) | LC-6A, LC-10AS | 228-11200-91 | kit | 25289 |
| Outlet Check Valve | LC-600, LC-9A, LC-10AD, LC-10AT | 228-18522-92, 228-32531-92 | ea. | 25282 |
| Outlet Check Valve | LC-10ADvp, LC-10ATvp | 228-34976-91 | ea. | 24983 |
| Outlet Check Valve | LC-20AD/AB XR | 228-45705-91, 228-45563-91 | ea. | 26427 |
| Sapphire Plunger | LC-10AS, LC-7A | 228-17019-93 | ea. | 25291 |
| Sapphire Plunger | LC-600, LC-9A, LC-10AD | 228-18523-91 | ea. | 25294 |
| Plunger Assembly, Ceramic | LC-10ADvp | 228-35601-91 | ea. | 25472 |
| Plunger Assembly, Ceramic | LC-10ATvp | 228-35009-92 | ea. | 25473 |
| Plunger Assembly, Sapphire | LC-10ADvp | 228-35601-92 | ea. | 24982 |

*Graphite-filled PTFE; best for organic solvents.

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion; best for buffers.



For even more options, visit www.restek.com/lcacc

