



Preventive Maintenance (PM) Kits for Waters HPLC Systems

Description	Model #	Similar to Waters part #	qty.	cat.#
2690/2695 Pump & Autosampler PM Kit Includes: sapphire plungers (2), seal wash plunger seals (4), head plunger seals (4), wash tube seals (4), battery, sparge diffuser, filter insert, face seals (4), solvent reservoir 10 µm filters (4), 250 µL syringe, check valve cartridges (2), wash tube seal, seal wash tube, PTFE washer, filter retainer, lower wash seal frit, needle wash frit, TFE washer, needle assembly, gold injector seals (2), stainless steel ferrule, compression screw	Alliance 2690, 2695	WAT270944	kit	25143
717 Autosampler PM Kit Includes: seal pack assembly, tube assembly (0.020" mm ID), needle, needle compression screw, 0.062 stainless steel ferrule, pre-column filter assembly, filter insert, 250 µL WISP syringe	717 Autosampler	WAT052669	kit	25145
600 Pump PM Kit Includes: PerformancePLUS cartridges (4), sparge diffusers (4), Super Seals (2), solvent reservoir 10 µm filters (4), sapphire plungers (2), reference valve button, valve disk spacer, valve disk, TFE ball plug, TFE seat, ruby ball, inlet tube body assembly manifold insert, insert seal, bellville washers (2), flat washer	600 Pump	WAT052675	kit	25144
616 Pump PM Kit Includes: sapphire plungers (2), sparge diffusers (4), solvent reservoir 10 µm filters (4), plunger seals (2), plunger wash seals (2), O-rings (2), backup rings (2), check valve cartridges (4)	616 Pump	WAT052672	kit	25146
515 Pump PM Kit Includes: PerformancePLUS check valves (4), sparge diffuser, solvent reservoir 10 µm filter, sapphire plungers (2), plunger seals (2), pivot inserts (2), pivot guides (2), washer (2), plunger springs (2), retaining rings (2)	515 Pump	WAT052587	kit	26519
1525 Pump PM Kit Includes: sapphire plungers (4), check valve cartridges (8), plunger seals (4), solvent reservoir 10 µm filters (2), reference valve button, valve disk spacer, valve disk	1525 Pump	201000114	kit	26430

Simplify Routine LC Maintenance Preventive Maintenance Kits

- Significant savings over instrument manufacturer prices.
- High-quality components in every kit.
- Wide range of options for HPLC systems and pumps.

www.restek.com/LC-Maintenance






Injector & Autosampler Supplies for Waters HPLC & ACQUITY UPLC® Systems

Description	Model #	Similar to Waters part #	qty.	cat. #
Needle Assembly, 0.010 mm ID	ACQUITY	700002644	ea.	25942
Sample Loop, 10 µL	ACQUITY, nanoACQUITY	430001326	ea.	25943
Wash Station O-Ring	ACQUITY, nanoACQUITY	700002572	ea.	25944
Metering Syringe, 100 µL	ACQUITY, nanoACQUITY	700002570	ea.	25945
Tube Assembly, Inject Outlet	ACQUITY	430001084	ea.	25946
Tube Assembly, System Outlet	ACQUITY	430001486	ea.	25947
Seal Pack Rebuild Kit, Without Seal Wash Tube Assembly Includes: wash tube seal, seal wash tube, PTFE washer, filter retainer, lower wash seal frit, needle wash frit, TFE washer, needle, injector seal (2), compression screw, ferrule	2690, 2695	WAT271019	kit	25495
Seal Pack With Needle Includes: seal pack assembly, pre-column filter assembly, needle, 0.020" ID tube assembly, compression screw, ferrule	717, LC Module 1	WAT045559	ea.	25496
Vespel Rotor Seal	1090, 7000, 7010, 7040, 7067	7010-039	ea.	25279
Vespel Rotor Seal	7125, 7126, 7725, 7725i, 9725	7125-047	ea.	25280
Isolation Seal	7010	7010-015	ea.	25281

Instrument Replacement Parts for Waters HPLC & ACQUITY UPLC® Systems



Pump Supplies for Waters HPLC & ACQUITY UPLC® Systems (Also see pages 310–311.)

Description	Model #	Similar to Waters part #	qty.	cat.#
Head Plunger Seal Kit Includes: head plunger seals (2), back-up washers (2)	ACQUITY, nanoACQUITY	700002599	kit	26428
Wash Seal	ACQUITY, nanoACQUITY	700002598	2-pk.	26429
Pump O-Ring	ACQUITY, nanoACQUITY	WAT076152	ea.	25951
In-Line Filter Assembly, Stainless Steel Frit	ACQUITY, nanoACQUITY	289003547	ea.	25952
Filter Frit Cartridge, Stainless Steel	ACQUITY	700002913	ea.	25953
Solvent Bottle Filters	ACQUITY	700003616	7-pk.	25954
Tube Assembly, SSV to In-line Filter	nanoACQUITY	430001470	ea.	25948
Tube Assembly, Transducer to Check Valve, BSM	ACQUITY, nanoACQUITY	430001121	ea.	25949
Tube Assembly, Transducer to Check Valve, QSM	ACQUITY	430002357	ea.	25950
Tube Assembly, SSV to I2V	ACQUITY	430001443	ea.	25941
Super Seal for Analytical Heads	M6KA, 1515, 1525, 510, 515, 590, 600, 610, LC Module 1	WAT022946	ea.	25374
Plunger Seal, Gold**	1515, 1525, 510, 515, 600, 610, LC Module 1	—	ea.	25375
Plunger Seal, Black*	1515, 1525, 510, 515, 600, 610, LC Module 1	WAT026613	ea.	25378
Plunger Seal, Gold for EF Heads**	510, 515, 590, 600, 610, LC Module 1	700002282	ea.	25380
Plunger Seal, Black for EF Heads*	510, 515, 590, 600, 610, LC Module 1	WAT026644	ea.	25379
Pump Seal Kit Includes: plunger seals (2), wash seals (2), O-rings (2), back-up rings (2)	616, 1525 Micro	WAT034515	kit	25497
Seal Wash Plunger Seals	2690, 2695, 2790, 2795	WAT270160, WAT271018	2-pk.	25386
Head Plunger Seals, Gold	2690, 2695, 2790, 2795	WAT270789, WAT270938	2-pk.	25387
Head Plunger Seals, Black	2690, 2695, 2790, 2795	WAT271066	2-pk.	25388
Seal Wash Face Seal	2690, 2695, 2790, 2795	WAT271017	ea.	25428
Wash Tube Seal Kit	2690, 2695, 2790, 2795	WAT270940	4-pk.	25429
Face Seals Replacement Kit	2690, 2695, 2790, 2795	WAT270163, WAT270939	4-pk.	26547



*Graphite-filled PTFE; best for organic solvents.

**Ultra-high molecular weight polyethylene (UHMWPE); increased resistance to abrasion; best for buffers.



Pump Supplies for Waters HPLC & ACQUITY UPLC® Systems (cont.)

Description	Model #	Similar to Waters part #	qty.	cat.#
Primary Check Valve	ACQUITY, nanoACQUITY	700002596	2-pk.	25955
Accumulator Check Valve (Double Ball & Seat)	ACQUITY, nanoACQUITY	700002968	2-pk.	25956
I2 Check Valve Cartridge	ACQUITY	700005165	ea.	25957
Inlet Check Valve Assembly	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT033679, WAT025214	ea.	25360
Inlet Check Valve Housing	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT025203	ea.	25361
Inlet Check Valve Rebuild Kit Includes: blue washer (2), retainer gasket (2), Kel-F insert (2), ruby ball (2), TFE seat (2), TFE washer (2)	M6KA, 510, 515, 590, 600, 610, LC Module 1	WAT060495	kit	25362
Outlet Check Valve Assembly (Actuator Style)	M6KA, 510, 590	WAT025028	ea.	25363
Outlet Check Valve Housing (Actuator Style)	M6KA, 510, 590	WAT025212	ea.	25364
Outlet Check Valve Rebuild Kit (Actuator Style) Includes: 2 µm cup filter (2), compression spring (2), actuator (2), TFE washer actuator seat (2), TFE washer (2)	M6KA, 510, 590	WAT026016	kit	25365
Outlet Check Valve Assembly (Ball & Seat Style)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT025216	ea.	25366
Outlet Check Valve Housing (Ball & Seat Style)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT025207	ea.	25367
Outlet Check Valve Rebuild Kit (Ball & Seat Style) Includes: 2 µm cup filter (2), ruby ball (2), ball seat (2), Kel-F insert (2), blue washer (2), washer (2)	1525EF, 510, 515, 600, 610, LC Module 1	WAT026014	kit	25368
Inlet Check Valve Assembly, 225 µL (Extended Flow)	1525EF, 510, 515, 600, 610, LC Module 1	WAT060307	ea.	25369
Check Valve Rebuild Kit (Extended Flow) Includes: 3/16" EF ruby ball (2), EF ball seat (2), EF guide seal (2), EF ball guide (2), EF housing seal (2)	510, 590, 600, LC Module 1	WAT088223	kit	25371
PerformancePLUS Check Valve Housing	1515, 1525, 510, 515, 525, 600, 610	—	ea.	25372
PerformancePLUS Check Valve Cartridge	1515, 1525, 2690, 2695, 2795, 510, 515, 600, 610	700000254	2-pk.	25370
Check Valve Cartridges	2690, 2695, 2790, 2795	WAT270941	2-pk.	25373
Check Valve Cartridge, Stainless Steel	616	WAT024960	ea.	25498

Instrument Replacement Parts for Waters HPLC & ACQUITY UPLC® Systems



Pump Supplies for Waters HPLC & ACQUITY UPLC® Systems (cont.)

Description	Model #	Similar to Waters part #	qty.	cat.#
Sapphire Plungers	ACQUITY, nanoACQUITY	700002600	2-pk.	25958
Sapphire Plunger	M6KA, 510, 590, 600, 610, LC Module 1	WAT025656, WAT069511	ea.	25381
Sapphire Plunger (Extended Flow)	1525EF, 510, 515, 590, 600, 610, LC Module 1	WAT060304	ea.	25382
Sapphire Plunger Assembly Kit Includes: sapphire plunger, pivot insert, pivot guide, washer, spring, retaining ring	515, 1515, 1525	WAS207069	kit	25384
Sapphire Plunger	616, 625, 626, 1525 Micro	WAT031788, WAT0270488	ea.	25420
Sapphire Plunger	2690, 2695, 2790, 2795, 2796	WAT270488, WAT2710678, WAT270959	ea.	25385
Draw-Off Tube Assembly	M6KA, 510, 515, 590, 600, 610	WAT060476	ea.	25392
Indicator Manifold Kit Includes: TFE ball plug, TFE seat, ruby ball, manifold insert, TFE insert seal, belleville washer (2), flat washer, inlet tube body assembly	M45, 501, 510, 515, 590, 600, 610	WAT060448	kit	25412
Reference Valve Rebuild Kit Includes: ref valve button, valve disk spacer, valve disk	600 Pump	WAT025746	kit	25492
Ferrule, Stainless Steel	—	WAT022330	ea.	25417
Compression Screws, Stainless Steel	—	WAT025313	2-pk.	25493
Battery	2690, 2695, 717, 486, 484	WAT080443	ea.	25494



25404



25405



25407



25408



25406



25410



25409



25423

25424



25427

Detector Supplies for Waters HPLC & ACQUITY UPLC® Systems

Description	Model #	Similar to Waters part #	qty.	cat.#
Xenon Lamp (without holder or mirror)	470	—	ea.	25404
Xenon Lamp	474	—	ea.	25405
Deuterium Lamp (UV/Vis)	484	WAT080357	ea.	25406
Deuterium Lamp (UV/Vis)	486	WAT080678	ea.	25407
Deuterium Lamp, Long Life (2,000 hours)*	486	WAT080678	ea.	25410
Deuterium Lamp	996, 2996	WAT052586	ea.	25408
Deuterium Lamp	2487, 2488	WAS081142	ea.	25409
Window Gasket	484, 486, 490	WAT080335	ea.	25423
Lamp Side Gasket	484, 486, 490	WAT080336	ea.	25424
Quartz Lens	486	WAT080687	ea.	25427

*Standard lamps have nominal 1,000-hour life.

More Choices • More Savings

Thomson SINGLE StEP™ Filter Vials

Save time and reduce waste—these eco-friendly vials replace syringe filtration and eliminate the need for separate vials, caps, syringes, and filters. Remove contaminants and particulates from your HPLC samples with just a single squeeze!

- Nylon, PES, PTFE & PVDF membranes in 0.2 µm and 0.45 µm porosities available.
- Fast and easy to use—only hand depression force is required.
- Compatible with both UHPLC and standard HPLC autosamplers.
- Pre-slit caps are color coded for easy identification of membrane and porosity.
- Low dead volume unit contains rugged polypropylene insert with 450 µL loading capacity.

Patent No.7,790,117



www.restek.com/singlestep