

Long Axis Rotating Mixer Shaker



Timing function




Six different modules

TECHNICAL DATA

Model	RML - 80	RML - 80Pro
Motor	Brush motor	Brushless motor
Control Mode	Knob	Digital display & Button
Type	Non-stop run	Non-stop run / Timing
Speed	20 - 80rpm	
Timing	/ 1min - 23h59min	
Power	25W	
Dimension (mm)	425*100*125	

MODULES

	Composite module	Centrifuge tube : 50ml*4+15ml*10+2ml*8 Blood collection tube : 10ml*10+2/5ml*8
	RML-801	Centrifuge tube : 50ml*12
	RML-802	Centrifuge tube : 15ml*23 Blood collection tube : 10ml*23
	RML-803	Centrifuge tube : 1.5/2.0ml*28 Blood collection tube : 2.0/5.0ml*28
	RML-804	Centrifuge tube : 15ml*16 Blood collection tube : 10ml*16
	RML-805	Centrifuge tube : 50ml*12

Rotary Rotating Mixer Shaker



Timing function



Large handling capacity



Adjustable angle



Brushless motor with strong power

TECHNICAL DATA

Model	RMO - 80	RMO - 80Pro
Motor	Brush motor	Brushless motor
Control Mode	Knob	Digital display & Button
Type	Non-stop run	Non-stop run / Timing
Speed	20 - 80rpm	
Angle Range	40° - 80°	
SPECS	Centrifuge tube : 50ml*8 or 15ml*16 or 1.5ml/2ml*24 Blood collection tube : 10ml*16 or 2ml/5ml*24	
Timing	/	1min - 23h59min
Allowed Ambient Humidity Temp.	5 ~ 40 C 80%RH	
Voltage	AC 100 - 220V, 50/60Hz	
Power	20W	
Dimension (mm)	230*180*265mm	